


Suspicion of peanut allergy

Is it allergy? Risk for severe reactions?



- **Ara h 2 – ImmunoCAP® f423 rAra h 2**
 - Often associated with systemic and severe reactions.
 - Associated with cross-reactivity to tree nuts, e.g. almond and brazil nut.
 - A peanut storage protein.
- **Ara h 1 – ImmunoCAP® f422 rAra h 1**
 - Often associated with clinical reactions.
 - Associated with cross-reactivity to nuts and legume seeds, e.g. lentil and pea.
 - A peanut storage protein.
- **Ara h 3 – ImmunoCAP® f424 rAra h 3**
 - Often associated with clinical reactions.
 - Associated with cross-reactivity to lupin and soybean.
 - A peanut storage protein.
- **Ara h 9 – ImmunoCAP® f427 rAra h 9**
 - Often associated with systemic and more severe reactions in addition to Oral Allergy Syndrome (OAS).
In peanut however, severe reactions to LTP are not well documented.
 - Stable to heat and digestion, risk for reactions also to cooked food.
 - A LTP (Lipid Transfer Protein).
- **Ara h 8 – ImmunoCAP® f352 rAra h 8**
 - Often associated with local reactions such as OAS.
 - A marker for birch-pollen related cross-reactions to peanut.
 - A heat labile protein, cooked food is often tolerated.
 - A PR-10 protein.
- **Birch (t3)**

Clinical presentation	Reason for requesting ImmunoCAP Allergen Components	Results of ImmunoCAP Allergens (kU/l)																																																
<ul style="list-style-type: none"> 17 year-old boy with increased problems of eczema, possible asthma. Consultation for a food allergy investigation. Eczema since 1 year of age, rhinitis during spring, occasions with breathing problems during past three years. Allergy to cat, grass and birch. Worsening of eczema the last six months. Family history. No known allergies in the family. SPT results Positive +3 for birch, grass and cat Positive +1 for dog, peanut, hazelnut and soy Negative for grass, wheat and soy. ImmunoCAP Specific IgE (kU/l) <table border="1"> <tr><td>Hazelnut</td><td>> 100</td></tr> <tr><td>Peanut</td><td>5.1</td></tr> <tr><td>Soy</td><td>2.0</td></tr> <tr><td>Rye</td><td>1.2</td></tr> <tr><td>Wheat</td><td>1.1</td></tr> <tr><td>Dog</td><td>1.9</td></tr> </table> 	Hazelnut	> 100	Peanut	5.1	Soy	2.0	Rye	1.2	Wheat	1.1	Dog	1.9	<p>Decided to perform a risk assessment of the peanut and hazelnut allergy.</p> <p>Test also for sensitization to soy in combination with a test for birch.</p>	<table border="1"> <tr><td>f13</td><td>Peanut</td><td>2.7</td></tr> <tr><td>f422</td><td>Ara h 1</td><td><0.1</td></tr> <tr><td>f423</td><td>Ara h 2</td><td><0.1</td></tr> <tr><td>f424</td><td>Ara h 3</td><td>< 0.1</td></tr> <tr><td>f352</td><td>Ara h 8</td><td>2.1</td></tr> <tr><td>f427</td><td>Ara h 9</td><td>< 0.1</td></tr> <tr><td>f17</td><td>Hazelnut</td><td>166</td></tr> <tr><td>f428</td><td>Cor a 1</td><td>80.6</td></tr> <tr><td>f425</td><td>Cor a 8</td><td>< 0.1</td></tr> <tr><td>f14</td><td>Soy</td><td>0.42</td></tr> <tr><td>f353</td><td>Gly m 4</td><td>38.5</td></tr> <tr><td>t3</td><td>Birch</td><td>206</td></tr> </table>	f13	Peanut	2.7	f422	Ara h 1	<0.1	f423	Ara h 2	<0.1	f424	Ara h 3	< 0.1	f352	Ara h 8	2.1	f427	Ara h 9	< 0.1	f17	Hazelnut	166	f428	Cor a 1	80.6	f425	Cor a 8	< 0.1	f14	Soy	0.42	f353	Gly m 4	38.5	t3	Birch	206
Hazelnut	> 100																																																	
Peanut	5.1																																																	
Soy	2.0																																																	
Rye	1.2																																																	
Wheat	1.1																																																	
Dog	1.9																																																	
f13	Peanut	2.7																																																
f422	Ara h 1	<0.1																																																
f423	Ara h 2	<0.1																																																
f424	Ara h 3	< 0.1																																																
f352	Ara h 8	2.1																																																
f427	Ara h 9	< 0.1																																																
f17	Hazelnut	166																																																
f428	Cor a 1	80.6																																																
f425	Cor a 8	< 0.1																																																
f14	Soy	0.42																																																
f353	Gly m 4	38.5																																																
t3	Birch	206																																																
<h3>Interpretation</h3> <ul style="list-style-type: none"> Indication of birch pollen induced hazelnut allergy as increased level of specific IgE to Cor a 1. Pollen-driven food allergy is in most cases associated with mild oral allergy syndrome (OAS). Special caution in this case: The very high IgE result of f17 Hazelnut indicates that other proteins might be involved in the sensitisation of this patient, risk for severe reactions should not be neglected Indication of birch pollen induced peanut allergy as increased levels of specific IgE to Ara h 8. Pollen-driven food allergy is in most cases associated with milder forms of OAS. Increased levels of specific IgE for Gly m 4 indicates a risk of severe reactions caused by soy-intake, particularly in high amounts and in combination with exercise during pollen seasons. 																																																		

Clinical presentation	Reason for requesting ImmunoCAP Allergen Components	Results of ImmunoCAP Allergens (kU/l)																				
<ul style="list-style-type: none"> 16 year-old girl on follow-up visit after a reaction caused by inhaled peanut. Swelling of the throat after inhaled peanut in an aeroplane. Eczema in early childhood, allergic to peanut since early childhood (local reaction in the mouth) and to furred animals since the age of 7 (rhinoconjunctivitis). Family history. Younger brother has a history of allergy to egg and milk, now tree-pollen allergic, as is the father. ImmunoCAP Specific IgE results: <table border="1"> <tr><td>Peanut</td><td>1.0 kU/l</td></tr> </table> 	Peanut	1.0 kU/l	<p>Decided to perform a risk assessment of the peanut allergy before deciding on a provocation test.</p>	<table border="1"> <tr><td>f13</td><td>Peanut</td><td>1.6</td></tr> <tr><td>f422</td><td>Ara h 1</td><td><0.1</td></tr> <tr><td>f423</td><td>Ara h 2</td><td>1.4</td></tr> <tr><td>f424</td><td>Ara h 3</td><td>< 0.1</td></tr> <tr><td>f352</td><td>Ara h 8</td><td>< 0.1</td></tr> <tr><td>f427</td><td>Ara h 9</td><td>< 0.1</td></tr> </table>	f13	Peanut	1.6	f422	Ara h 1	<0.1	f423	Ara h 2	1.4	f424	Ara h 3	< 0.1	f352	Ara h 8	< 0.1	f427	Ara h 9	< 0.1
Peanut	1.0 kU/l																					
f13	Peanut	1.6																				
f422	Ara h 1	<0.1																				
f423	Ara h 2	1.4																				
f424	Ara h 3	< 0.1																				
f352	Ara h 8	< 0.1																				
f427	Ara h 9	< 0.1																				
<h3>Interpretation</h3> <ul style="list-style-type: none"> The sensitisation to the peanut allergens is induced by the Ara h 2 protein, a storage protein often associated with severe reactions. Even though the peanut specific IgE level is low in this case there is an increased risk of severe reactions. 																						